TURB - ASO Anti-Estreptolisina O

Finalidade:

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1. Stevens DL. Invasive group A streptococcal disease. Clin Infect Dis 1992; 14: 2-11.
2. Immunology and Serology in Laboratory Medicine, 2nd edition. Turgeon mL. Mosby, 1996.
3. Arquivos da Ebram.

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	Calibrator (Slope AVG)
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	0
Unit:	UI/mL
ANALYSIS	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name: H2O
Post Dil Factor: 2	Conc. Factor: NO
Sample:	Cycle: 1
Volume:	Dil: 20ul
Reagent:	Cycle: 1
Volume:	240ul
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil: 2
CALCULATION	
Sample Limit:	No
Point:	
Reaction Direction:	Increase
Check:	On
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	0
High:	800
Norm Range Low:	0
High:	200
Number of Steps:	1
Calculation Step A:	ENDPOINT
Readings First:	1
Last:	8
Reaction Limit:	NO
Point:	
Calib. Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	NO
Calibrator Pos:	(*)
STD1:	STD2:
STD3:	STD4:
STD5:	STD6:
STD7:	STD8:
Calc. Model:	
Correction STD:	
Replicate:	DUPL
Deviation:	10%
CONTROL	
CS1 -	
CS2 -	
CS3 -	

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Test Name: ASO	Test: ASO		
Test Bar Code:			
Test Type: End Point	Curve Type: Blanked Linear		
Units: UI/ml	N° of Decimal Places: 0		
Primary Wavelength: 540	Secondary Wavelength:		
Read Time Interval: 20	Sample Blank: No		
Factor:			
Calibration Interval: 168 hours			
Normalization Interval:			
N° of Calibrations: 2	N° of Replicates: 2		
Low Blank A Limit: -0.100	High Blank A Limit: 2.000		
Low A Limit: -0.100	High A Limit: 2.500		
Low Normal: 0	High Normal: 2.00		
Linearity Limit: 800	Curve S.D Limit: 8.0		
Test Name: ASO	Test: ASO		
Test Bar Code:			
Sample Volume: 3ul	Sample Diluent:		
Reagent Dilution Ratio: 1	Dilution Ratio: 1		
Reagent Dilution:			
Reagent Volume	Bar Code	Diluent Volume	Lag Time
Reagent 1	300	AS14	120
Reagent 2			
Reagent 3			
Reagent 4			
Control:			

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	Calibrator (Slope AVG)
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	0
Unit:	UI/mL
ANALYSIS	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name: H2O
Post Dil Factor: 2	Conc. Factor: NO
Sample:	Cycle: 1
Volume:	Dil: 20ul
Reagent:	Cycle: 1
Volume:	240ul
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil: 2
CALCULATION	
Sample Limit:	No
Point:	
Reaction Direction:	Increase
Check:	On
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	0
High:	800
Norm Range Low:	0
High:	200
Number of Steps:	1
Calculation Step A:	ENDPOINT
Readings First:	1
Last:	8
Reaction Limit:	NO
Point:	
Calib. Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	NO
Calibrator Pos:	(*)
STD1:	STD2:
STD3:	STD4:
STD5:	STD6:
STD7:	STD8:
Calc. Model:	
Correction STD:	
Replicate:	
Deviation:	
CS1 -	
CS2 -	
CONTROL	
CS3 -	

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Name:	ASO
Abbr Name:	ASO
Mode:	TwoPoint
Wavelength:	546 nm
Units:	UI/mL
Decimals:	0
Low Conc:	0
High Conc:	800
Calibrator Name:	Pad ASO
Repeat:	2
Number:	1
Concentration:	(*)
Interval:	(*)
Cut off:	No
Prozone Check:	No
Ref Male Low:	0
Ref Male High:	200
Ref Female Low:	0
Ref Female High:	200
Ref Ped Low:	0
Ref Ped High:	150
Control 1:	
Control 2:	
Control 3:	
Correlat Factor:	1.000
Correlat Offset:	0.000
DUAL MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
SAMPLE	
Normal Volume:	
Reagent Volume:	
R2 Bottle:	
Normal Volume:	
Reagent Volume:	
Predilution:	
Slope Blank:	
Delay Min Time:	
Linearity Limit:	
Point One Two:	
Incubation Time:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
R ABS Deviation:	
Reagent Blank:	
Cal Low Limit:	
Cal High Limit:	
Factor:	
MONO MODE	
Name:	ASO
Sample Blank:	NO
R1 Bottle:	25mL
Normal Volume:	300ul
Reagent Volume:	300ul
SAMPLE	
Normal Volume:	3ul
Reagent Volume:	2ul
Delay Min Time:	
Linearity Limit:	
Predilution:	NO
Point One Two:	12.110
Incubation Time:	
Low Absorbance:	-0.100
High Absorbance:	3.000
R ABS L Limit:	-0.100
R ABS H Limit:	1.000
R ABS Deviation:	2.000
Reagent Blank:	YES (**)
Cal Low Limit:	(*)
Cal High Limit:	(*)
Factor:	(*)

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