TURB - PCR Proteína C Reativa

Finalidade:

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1. Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER, WB Saunders Co., 1999.
2. Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
3. Arquivos Ebram

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	Calibrator (Slope AVG)
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	0
Unit:	mg/L
ANALYSIS	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name:
Post Dil Factor: 2	Conc. Factor: NO
Sample:	Cycle: 1
Volume: 3 uL	Dil: 10 uL
Reagent:	Cycle: 1
Volume: 240 uL	
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil:
CALCULATION	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	0
High:	150
Norm Range Low:	0
High:	5
Number of Stans:	1
Calculation Stan A:	ENDPOINT
Readings First:	1
Last:	8
Reaction Limit:	NO
Point:	
Calib. Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	
Calibrator Pos:	(*)
STD1:	STD2:
STD3:	STD4:
STD5:	STD6:
STD7:	STD8:
Calc. Model:	
Correction STD:	
Replicate:	DUPL
Deviation:	10%
CONTROL	
CS1 -	
CS2 -	
CS3 -	

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Test Name: PCR	Test: PCR		
Test Bar Code:	Curve Type: Blanked Linear		
Test Type: Endpoint	N° of Decimal Places: 1		
Units: mg/L	Primary Wavelength: 540		
Secondary Wavelength:	Read Time Interval: 20		
Sample Blank: No	Factor:		
Calibration Interval:	Normalization Interval:		
N° of Calibrations: 2	N° of Replicates: 2		
Low Blank A Limit: -0.100	High Blank A Limit: 2.500		
Low A Limit: -0.100	High A Limit: 2.000		
Low Normal: 0.0	High Normal: 5.0		
Linearity Limit: 150	Curve S.D Limit: 15.0		
Test Name: PCR	Test: PCR		
Test Bar Code:			
Sample Volume: 3 uL	Sample Diluent:		
Reagent Dilution Ratio: 1	Dilution Ratio: 1		
Reagent Dilution:			
Reagent Volume	Bar Code	Diluent Volume	Lag Time
Reagent 1	300 uL	PCR1A	120
Reagent 2			
Reagent 3			
Reagent 4			
Control:			

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	Calibrator (Slope AVG)
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	0
Unit:	mg/L
ANALYSIS	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name:
Post Dil Factor: 2	Conc. Factor: NO
Sample:	Cycle: 1
Volume: 3 uL	Dil: 10 uL
Reagent:	Cycle: 1
Volume: 240 uL	
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil:
CALCULATION	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	0
High:	150
Norm Range Low:	0
High:	5
Number of Stans:	1
Calculation Stan A:	ENDPOINT
Readings First:	1
Last:	8
Reaction Limit:	NO
Point:	
Calib. Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	
Calibrator Pos:	(*)
STD1:	STD2:
STD3:	STD4:
STD5:	STD6:
STD7:	STD8:
Calc. Model:	
Correction STD:	
Replicate:	
Deviation:	
CS1 -	
CS2 -	
CONTROL	
CS3 -	

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Name:	PCR
Abbr Name:	
Mode:	TwoPoint
Wavelength:	546
Units:	mg/L
Decimals:	1
Low Conc:	0
High Conc:	150
Calibrator Name:	(*)
Repeat:	2
Number:	1
Concentration:	(**)
Interval:	(*)
Cut off:	No
Prozone Check:	No
Ref Male Low:	0
Ref Male High:	5
Ref Female Low:	0
Ref Female High:	5
Ref Ped Low:	0
Ref Ped High:	5
Control 1:	
Control 2:	
Control 3:	
Correlat Factor:	1.000
Correlat Offset:	0.000
DUAL MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
SAMPLE	
Normal Volume:	
Reagent Volume:	
R2 Bottle:	
Normal Volume:	
Reagent Volume:	
Predilution:	
Slope Blank:	
Delay Min Time:	
Linearity Limit:	
Point One Two:	
Incubation Time:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
R ABS Deviation:	
Reagent Blank:	
Cal Low Limit:	
Cal High Limit:	
Factor:	
MONO MODE	
Name:	PCR
Sample Blank:	No
R1 Bottle:	25ml
Normal Volume:	300ul
Reagent Volume:	300ul
SAMPLE	
Normal Volume:	2ul
Reagent Volume:	1ul
Delay Min Time:	
Linearity Limit:	
Predilution:	No
Point One Two:	12.110
Incubation Time:	
Low Absorbance:	-0.100
High Absorbance:	3.000
R ABS L Limit:	-0.100
R ABS H Limit:	1.000
R ABS Deviation:	2.000
Reagent Blank:	YES (***)
Cal Low Limit:	(*)
Cal High Limit:	(*)
Factor:	(***)

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