TURB FR - Fator Reumatóide FR

Finalidade:

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	D-R-S-SR1
Calibration Mode:	STD NONLINEAR (ou LOGIT/LOG4)
Reagent Blank:	NO BLANK
Cleaner:	NO
Wavelength:	800nm
Decimal Position:	1
Unit:	U/ml
ANALYSIS	
Dilution Name:	H2O
Factor:	NO
Time:	NO
Post Dil Factor:	2.00
Sample:	Cycle: 2
Volume:	Dil: 0.0
Reagent:	Cycle: 1
Volume:	Dil: 0.0
Start Reagent 1:	Cycle: 1
Volume:	Dil: 0.0
Start Reagent 2:	Cycle: 1
Volume:	Dil: 0.0
CALCULATION	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.000
Offset:	0.000
Test Range Low:	ON
High:	ON
Norm Range Low:	0
High:	20
Number of Steps:	1
Calculation Stan A:	ENDPOINT
Readings First:	2
Last:	30
Reaction Limit:	
Point:	
Calib Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low:	NO
High:	NO
Blank Range Low:	NO
High:	NO
Factor:	
Calibrator Pos:	(**)
STD1:	(*)
STD2:	(*)/2
STD3:	(*)/4
STD4:	(*)/8
STD5:	(*)/16
STD6:	0.0
STD7:	
STD8:	
Calc Model:	LOGIT/LOG4
Correction STD:	NO
Replicate:	SINGLE
Deviation:	10%
CONTROL	
CS1 - NO	
CS2 - NO	
CS3 - NO	

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Test Name:	Fator Reumatóide	Test: FR		
Test Bar Code:				
Test Type:	Endpoint	Curve Type: 2-Inv5		
Units:	U/ml	Nº of Decimal Places:		
Primary Wavelength:	600nm	Secondary Wavelength:		
Read Time Interval:	20	Sample Blank:		
Factor:	1.0000			
Calibration Interval:	999			
Normalization Interval:	999			
Nº of Calibrators:	6	Nº of Replicates:		
Low Blank A Limit:	-0.100	High Blank A Limit:		
Low A Limit:	-0.1000	High A Limit:		
Low Normal:	0	High Normal:		
Linearity Limit:		Curve S.D. Limit:		
Test Name:	Fator Reumatóide	Test: FR		
Test Bar Code:				
Sample Volume:	5ul	Sample Diluent:		
Reagent Dilution Ratio:	1	Predilution Ratio:		
Reagent Dilution:				
	Reagent Volume	Bar Code	Diluent Volume	Lag Time
Reagent 1	170	FR1		30sec
Reagent 2	100	FR2		450sec
Reagent 3				
Reagent 4				
Control:				

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	D-R-S-SR1
Calibration Mode:	LOGIT/LOG4
Reagent Blank:	NO BLANK
Cleaner:	NO
Wavelength:	800nm
Decimal Position:	1
Unit:	U/ml
ANALYSIS	
Dilution Name:	H2O
Factor:	NO
Time:	NO
Post Dil Factor:	2.00
Sample:	Cycle: 2
Volume:	Dil: 0.0
Reagent:	Cycle: 1
Volume:	Dil: 0.0
Start Reagent 1:	Cycle: 1
Volume:	Dil: 0.0
Start Reagent 2:	Cycle: 1
Volume:	Dil: 0.0
CALCULATION	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.000
Offset:	0.000
Test Range Low:	ON
High:	ON
Norm Range Low:	0
High:	20
Number of Steps:	1
Calculation Stan A:	ENDPOINT
Readings First:	2
Last:	30
Reaction Limit:	
Point:	
Calib Interval:	
CALIBRATION	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low:	NO
High:	NO
Blank Range Low:	NO
High:	NO
Factor:	
Calibrator Pos:	(**)
STD1:	(**)
STD2:	(**)
STD3:	(**)
STD4:	(**)
STD5:	(**)
STD6:	(**)
STD7:	(**)
STD8:	(**)
Calc Model:	
Correction STD:	NO
Replicate:	SINGLE
Deviation:	10%
CONTROL	
CS1 - NO	
CS2 - NO	
CS3 - NO	

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Name:	FR
Abbr Name:	FR
Mode:	Endpoint
Wavelength:	620nm
Units:	U/ml
Decimals:	1
Low Conc:	0.0 U/ml
High Conc:	120.0 U/ml
Calibrator Name:	Calib FR
Renal:	1
Number:	6
Concentration:	(*)
Interval:	(*)
Cut off:	NO
Prozone Check:	NO
Ref Male Low:	(*)
Ref Male High:	(*)
Ref Female Low:	(*)
Ref Female High:	(*)
Ref Pad Low:	(*)
Ref Pad High:	(*)
Control 1:	
Control 2:	
Control 3:	
Correlat Factor:	1.000
Correlat Offset:	0.000
DUAL MODE	
Name:	FR
Sample Blank:	NO
R1 Bottle:	25ml
Normal Volume:	220ul
Runn Volume:	224ul
SAMPLE	
Normal Volume:	8ul
Runn Volume:	4ul
R2 Bottle:	5ml
Normal Volume:	130ul
Runn Volume:	130ul
Predilution:	NO
Slope Blank:	
Delay Min Time:	
Linearity Limit:	
Point One, Two:	
Incubation Time:	6.5 min
Low Absorbance:	-0.100
High Absorbance:	3.000
R ABS L Limit:	-0.100
R ABS H Limit:	1.000
R ABS Deviation:	
Reagent Blank:	YES (**)
Cal Low Limit:	(*)
Cal High Limit:	(*)
Factor:	(**)
MONO MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Runn Volume:	
SAMPLE	
Normal Volume:	
Runn Volume:	
Delay Min Time:	
Linearity Limit:	
Predilution:	
Point One, Two:	
Incubation Time:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
R ABS Deviation:	
Reagent Blank:	
Cal Low Limit:	
Cal High Limit:	
Factor:	

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