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## Finalidade:

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1. Herbeth B, Siest G, Henny J. High-sensitivity C- reactive protein (CRP) reference intervals in the elderly. Clin Chem Lab Med 2011; 39: 1169-70
2. Price CP, Trull AK, Berry D, Gorman EG. Development and validation of a particle-enhanced turbidimetric immunoassay for C-reactive protein. J Immunol Methods 1987; 99: 205-211
3. Friedman and Young. Effects of disease on clinical laboratory tests. 3th ed. AACCPress, 1997.
4. Pearson TA, Mensah GA, Alexander RW, et al. Circulation, 2003, 107:499-511.
4. Arquivos Ebram

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<b>GENERAL</b>	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	SLOPE/AVG ou CALIBRATOR
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	2
Unit:	mg/L
<b>ANALYSIS</b>	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name: H2O
Post Dil Factor: 2.00	Conc. Factor: NO
Sample:	Cycle: 1
Volume:	Dil: 10 uL
Reagent:	Cycle: 1
Volume:	
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil:
<b>CALCULATION</b>	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.000
Offset:	0.000
Test Range Low:	0
High:	15
Norm Range Low:	NO
High:	NO
Number of Stans:	1
Calculation Step A:	ENDPOINT
Readline First:	T1
Last:	15
Reaction Limit:	NO
Point:	
Calib. Interval:	
<b>CALIBRATION</b>	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	
Calibrator Pos:	(**)
STD1:	(*)
STD2:	
STD3:	
STD4:	
STD5:	
STD6:	
STD7:	
STD8:	
Calc. Model:	
Correction STD:	NO
Replicate:	DUPL
Deviation:	10%
<b>CONTROL</b>	
CS1 - NO	
CS2 - NO	
CS3 - NO	

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Test Name: PCR-US	Test: PCR-US		
Test Bar Code:			
Test Type: Endpoint	Curve Type: Blanket Linear		
Units: mg/L	N° of Decimal Places: 2		
Primary Wavelength: 540	Secondary Wavelength: No		
Read Time Interval: 20	Sample Blank: No		
Factor:			
Calibration Interval: 999 horas			
Normalization Interval: 999 horas			
N° of Calibrators: 2	N° of Replicates: 2		
Low Blank A Limit: -0.100	High Blank A Limit: 2.500		
Low A Limit: -0.100	High A Limit: 2.000		
Low Normal: 0.00	High Normal: 8.50		
Linearity Limit: 15.00	Curve S.D Limit: 15.00		
Test Name: PCR-US	Test: PCR-US		
Test Bar Code:			
Sample Volume: 5 uL	Sample Diluent:		
Reagent Dilution Ratio: 1	Dilution Ratio: 1		
Reagent Dilution:			
Reagent Volume	Bar Code	Diluent Volume	Lag Time
Reagent 1	300 uL	PCRU	300
Reagent 2			
Reagent 3			
Reagent 4			
Control:			

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<b>GENERAL</b>	
Measurement Mode:	ABSORB
Reaction Mode:	R-S
Calibration Mode:	SLOPE/AVG
Reagent Blank:	REAG/DIL
Cleaner:	NO
Wavelength:	550 nm
Decimal Position:	2
Unit:	mg/L
<b>ANALYSIS</b>	
Dilution Name:	STD:
Factor:	Main STD:
Time:	Sample Dil Name: H2O
Post Dil Factor: 2.00	Conc. Factor: NO
Pos:	
STD1:	STD2:
STD3:	STD4:
STD5:	STD6:
STD7:	STD8:
Sample:	Cycle: 1
Volume:	Dil: 10 uL
Reagent:	Cycle: 1
Volume:	
Start Reagent 1:	Cycle:
Volume:	Dil:
Start Reagent 2:	Cycle:
Volume:	Dil:
<b>CALCULATION</b>	
Sample Limit:	NO
Point:	
Reaction Direction:	INCREASE
Check:	ON
Conversion Factor:	1.000
Offset:	0.000
Test Range Low:	0
High:	15
Norm Range Low:	NO
High:	NO
Number of Stans:	1
Calculation Step A:	ENDPOINT
Readline First:	T1
Last:	15
Reaction Limit:	NO
Point:	
Calib. Interval:	
<b>CALIBRATION</b>	
Reagent Blank:	ON REQUEST
Reagent Range:	
Reagent Range Low: NO	Blank Range Low: NO
High: NO	High: NO
Factor:	
Calibrator Pos:	(**)
STD1:	(*)
STD2:	
STD3:	
STD4:	
STD5:	
STD6:	
STD7:	
STD8:	
Calc. Model:	
Correction STD:	NO
Replicate:	DUPL
Deviation:	10%
<b>CONTROL</b>	
CS1 - NO	
CS2 - NO	
CS3 - NO	

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Name:	PCR-US
Abbr Name:	
Mode:	TwoPoint
Wavelength:	546 nm
Units:	mg/L
Decimals:	2
Low Conc:	0
High Conc:	15
Calibrator Name:	PCR Ultra
Repeat:	2
Number:	1
Concentration:	(*)
Interval:	
Cut off:	NO
Prozone Check:	NO
Ref Male Low:	(*)
Ref Male High:	(*)
Ref Female Low:	(*)
Ref Female High:	(*)
Ref Ped Low:	(*)
Ref Ped High:	(*)
Control 1:	
Control 2:	
Control 3:	
Correlat Factor:	1.000
Correlat Offset:	0.000
<b>DUAL MODE</b>	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
<b>SAMPLE</b>	
Normal Volume:	
Reagent Volume:	
R2 Bottle:	
Normal Volume:	
Reagent Volume:	
<b>Predilution</b>	
Sample Blank:	
Delay Min Time:	
Linearity Limit:	
Point One, Two:	
Incubation Time:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
R ABS Deviation:	
Reagent Blank:	
Cal Low Limit:	
Cal High Limit:	
Factor:	
<b>MONO MODE</b>	
Name:	PCR-US
Sample Blank:	No
R1 Bottle:	25mL
Normal Volume:	225uL
Reagent Volume:	225uL
<b>SAMPLE</b>	
Normal Volume:	3uL
Reagent Volume:	2uL
Delay Min Time:	
Linearity Limit:	
Predilution:	
Point One, Two:	12.304 sec
Incubation Time:	
Low Absorbance:	-0.100
High Absorbance:	3.000
R ABS L Limit:	-0.100
R ABS H Limit:	1.000
R ABS Deviation:	2.000
Reagent Blank:	Yes
Cal Low Limit:	(*)
Cal High Limit:	(*)
Factor:	(**)

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