QUIMIHA1C HEMOGLOBINA GLICOSILADA

Finalidade:

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1. Goldstein, D.E. et al, Diabetes Care. 27(7):1761-73 (2004)
2. United Kingdom Prospective study, Lancet 352: 837-53 (1998)
3. Little, R. et al., Clin Chemistry, 47: 1985-1992 (2001)
4. American Diabetes Association. Clinical practice recommendation: standards of medical care for patients with diabetes mellitus. Diab Care 22 (supp): S32-41 (1999)
5. Goldstein et al, Clin Chem 32: B64-B70 (1986)
6. Hoelzel W et al. IFCC reference system for measurement of hemoglobin A1c in human blood and the national standardization schemes in the United States, Japan and Sweden: a method-comparison study. Clin Chem 2004;50: 166-74

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| | | |
|-------------------------|-----------------------|-------|
| GENERAL | | |
| Measurement Mode: | ABSORB | |
| Reaction Mode: | R-S-SR1-SR2 | |
| Calibration Mode: | Lin Rear | |
| Reagent Blank: | NO Blank | |
| Cleaner: | No | |
| Wavelength: | 700 nm | |
| Decimal Position: | 2 | |
| Unit: | % | |
| ANALYSIS | | |
| Dilution Name: | STD: | |
| Factor: | Main STD: | |
| Time: | Sample Dil Name: | |
| Post Dil Factor NO: | Conc. Factor NO: | |
| Sample: | Circle: | 1 |
| Volume: | Dil: | 0.0ul |
| Reagent: | Circle: | 1 |
| Volume: | | |
| Star Reagent 1: | Circle: | 3 |
| Volume: | Dil: | 0.0ul |
| Star Reagent 2: | Circle: | 13 |
| Volume: | Dil: | 0.0ul |
| CALCULATION | | |
| Sample Limit: | No | |
| Point: | | |
| Reaction Direction: | Increase | |
| Check: | On | |
| Conversion Factor: | 1.0000 | |
| Offset: | 0.0000 | |
| Test Range Low: | Off | |
| High: | Off | |
| Norm Range Low: | No | |
| High: | No | |
| Number of Steps: | 1 | |
| Calculation Step A: | Endpoint | |
| Readings First: | 12 | |
| Last: | 20 | |
| Reaction Limit: | | |
| Point: | | |
| Calib. Interval: | | |
| CALIBRATION | | |
| Reagent Blank: | On request | |
| Reagent Range: | | |
| Reagent Range Low: -0.1 | Blank Range Low: -0.1 | |
| High: 0.4 | High: 0.4 | |
| Factor: | | |
| Calibrator Pos: | * | |
| STD1: | STD2: | ** |
| STD3: | STD4: | No |
| STD5: | STD6: | No |
| STD7: | STD8: | No |
| Calc. Model: | | |
| Correction STD: | | |
| Replicate: | Dual | |
| Deviation: | No | |
| CONTROL | | |
| CS1 - *** | | |
| CS2 - *** | | |
| CS3 - No | | |

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| | | | | |
|-------------------------|-----------------------|----------|----------------|----------|
| Test Name: | Test: | | | |
| Test Bar Code: | | | | |
| Test Type: | Curve Type: | | | |
| Unit: | Nº of Decimal Places: | | | |
| Primary Wavelength: | Secondary Wavelength: | | | |
| Read Time Interval: | Sample Blank: | | | |
| Factor: | | | | |
| Calibration Interval: | | | | |
| Normalization Interval: | | | | |
| Nº of Calibrators: | Nº of Replicates: | | | |
| Low Blank A Limit: | High Blank A Limit: | | | |
| Low A Limit: | High A Limit: | | | |
| Low Normal: | High Normal: | | | |
| Linearity Limit: | Curve S.D Limit: | | | |
| Test Name: | Test: | | | |
| Test Bar Code: | | | | |
| Sample Volume: | Sample Diluent: | | | |
| Reagent Dilution Ratio: | Predilution Ratio: | | | |
| Reagent Dilution: | | | | |
| | Reagent Volume | Bar Code | Diluent Volume | Lag Time |
| Reagent 1 | | | | |
| Reagent 2 | | | | |
| Reagent 3 | | | | |
| Reagent 4 | | | | |
| Comments: | | | | |

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| | |
|--------------------|--|
| Name: | |
| Abbr. Name: | |
| Mode: | |
| Wavelength: | |
| Unit: | |
| Decimals: | |
| Low Conc: | |
| High Conc: | |
| Calibrator Name: | |
| Repeat: | |
| Number: | |
| Concentration: | |
| Interval: | |
| Cut off: | |
| Prozone Check: | |
| Ref Male Low: | |
| Ref Male High: | |
| Ref Female Low: | |
| Ref Female High: | |
| Ref Pad Low: | |
| Ref Pad High: | |
| Control 1: | |
| Control 2: | |
| Control 3: | |
| Correl Factor: | |
| Correl Offset: | |
| DUAL MODE | |
| Name: | |
| Sample Blank: | |
| R1 Bottle: | |
| Normal Volume: | |
| Run Volume: | |
| SAMPLE | |
| Normal Volume: | |
| Run Volume: | |
| R2 Bottle: | |
| Normal Volume: | |
| Run Volume: | |
| Predilution: | |
| Slope Blank: | |
| Delay Min Time: | |
| Linearity Limit: | |
| Point One Two: | |
| Incubation Time: | |
| Low Absorbance: | |
| High Absorbance: | |
| R ABS L Limit: | |
| R ABS H Limit: | |
| Substr Denatation: | |
| Reagent Blank: | |
| R ABS Deviation: | |
| Cal Low Limit: | |
| Cal High Limit: | |
| Factor: | |
| MONO MODE | |
| Name: | |
| Sample Blank: | |
| R1 Bottle: | |
| Normal Volume: | |
| Run Volume: | |
| SAMPLE | |
| Normal Volume: | |
| Run Volume: | |
| Delay Min Time: | |
| Linearity Limit: | |
| Predilution: | |
| Incubation Time: | |
| Point One Two: | |
| Low Absorbance: | |
| High Absorbance: | |
| R ABS L Limit: | |
| R ABS H Limit: | |
| Substr Denatation: | |
| Reagent Blank: | |
| R ABS Deviation: | |
| Cal Low Limit: | |
| Cal High Limit: | |
| Factor: | |
| (*) | |
| (**) | |
| (***) | |
| (#) | |

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| | |
|------------------------|--|
| Name: | |
| Imunoassav: | |
| Chemicalv Test: | |
| Reagent: | |
| Inverse Chemicalv: | |
| % Sample Volume: | |
| Wavelength: | |
| Bichromatic Chemicalv: | |
| Bichromatic Factor: | |
| R1: | |
| R2: | |
| Re Limit 1: | |
| Re Limit 2: | |
| Deviation Limit: | |
| Delay Time: | |
| Blank Type: | |
| Incubate: | |
| % Reagent Volume: | |
| 2nd Reagent: | |
| 2nd Reagent Volume: | |
| A2 Delay: | |
| Unit: | |
| Unit Factor: | |
| Decimal Point: | |
| RBI Low: | |
| RBI High: | |
| Range Low: | |
| Range High: | |
| Calibration Factor: | |
| Reagent Rate: | |
| Standard Value: | |
| Normal Low: | |
| Normal High: | |
| Slope: | |
| Intersect: | |
| C1*10E-6: | |
| C2*10E-6: | |
| D1*10E-6: | |
| Delta R: | |
| Linear Factor: | |
| Final Limit: | |
| Endpoint Limit: | |
| (*) | |
| (**) | |

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| | | |
|--------------------|---------------|----------------|
| Measure Unit: | Decimals: | Reaction Type: |
| Volume (µl): | Reagent 1: | Urine: |
| Abs Range (m Abs): | Min: | Reagent 2: |
| Linearity Limit: | | Max: |
| Contaminating: | | Reagent Blank: |
| Filter 1: | | Differential: |
| Time (sec): | Mix 1: | Filter 2: |
| | Mix 2: | |
| Measurement Type: | Incubation 1: | Lag Phase: |
| Normal Range: | Incubation 2: | Measure: |
| | | Factor: |
| Age (yrs) | Homens | Mulheres |
| Below 10: | Min. | Max. |
| From 10 to 80: | | |
| Over: | | |

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