QUIMIHA1C HEMOGLOBINA GLICOSILADA

Finalidade:

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1. Goldstein, D.E. et al, Diabetes Care. 27(7):1761-73 (2004)
2. United Kingdom Prospective study, Lancet 352: 837-53 (1998)
3. Little, R. et al., Clin Chemistry, 47: 1985-1992 (2001)
4. American Diabetes Association. Clinical practice recommendation: standards of medical care for patients with diabetes mellitus. Diab Care 22 (supp): S32-41 (1999)
5. Goldstein et al, Clin Chem 32: B64-B70 (1986)
6. Hoelzel W et al. IFCC reference system for measurement of hemoglobin A1c in human blood and the national standardization schemes in the United States, Japan and Sweden: a method-comparison study. Clin Chem 2004;50:166-74

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GENERAL	
Measurement Mode:	ABSORB
Reaction Mode:	R-S-SR1-SR2
Calibration Mode:	Lin Rear
Reagent Blank:	NO Blank
Cleaner:	No
Wavelength:	700 nm
Decimal Position:	2
Unit:	%
ANALYSIS	
Dilution Name:	STD
Factor:	Main STD
Time:	Sample Dil Name
Post Dil Factor NO:	Conc. Factor NO
Sample:	Circle: 1
Volume:	Dil: 0.0ul
Reagent:	Circle: 1
Volume:	112ul
Start Reagent 1:	Circle: 3
Volume:	Dil: 0.0ul
Start Reagent 2:	Circle: 13
Volume:	Dil: 0.0ul
CALCULATION	
Sample Limit:	No
Point:	
Reaction Direction:	Increase
Check:	On
Conversion Factor:	1.0000
Offset:	0.0000
Test Range Low:	Off
High:	Off
Norm Range Low:	No
High:	No
Number of Slats:	1
Calculation Stop A:	Endpoint
Readings First:	12
Last:	20
Reaction Limit:	
Point:	
Calib Interval:	On request
CALIBRATION	
Reagent Blank:	
Reagent Range:	
Reagent Range Low: -0.1	Blank Range Low: -0.1
High: 0.4	High: 0.4
Factor:	
Calibrator Pos:	*
STD1:	**
STD2:	**
STD3:	No
STD4:	No
STD5:	No
STD6:	No
STD7:	No
STD8:	No
Calc Model:	
Correction STD:	
Replicate:	Dual
Deviation:	No
CONTROL	
CS1 - **	
CS2 - **	
CS3 - No	

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Disponemos do protocolo que executa o Lise on board. Contate o SAC

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Test Name:	Test:
Test Bar Code:	
Test Type:	Curve Type:
Units:	Nº of Decimal Places:
Primary Wavelength:	Secondary Wavelength:
Read Time Interval:	Sample Blank:
Factor:	
Calibration Interval:	
Normalization Interval:	
Nº of Calibrators:	Nº of Replicates:
Low Blank A Limit:	High Blank A Limit:
Low A Limit:	High A Limit:
Low Normal:	High Normal:
Linearity Limit:	Curve S.D. Limit:
Test Name:	Test:
Test Bar Code:	
Sample Volume:	Sample Dilution:
Reagent Dilution Ratio:	Predilution Ratio:
Reagent Dilution:	
	Reagent Volume
	Bar Code
	Diluent Volume
	Lag Time
Reagent 1	
Reagent 2	
Reagent 3	
Reagent 4	
Controls:	

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Name:	
Abbr Name:	
Mode:	
Wavelength:	
Units:	
Decimals:	
Low Conc:	
High Conc:	
Calibrator Name:	
Repeat:	
Number:	
Concentration:	
Interval:	
Cut off:	
Prozone Check:	
Ref Male Low:	
Ref Male High:	
Ref Female Low:	
Ref Female High:	
Ref Pad Low:	
Ref Pad High:	
Control 1:	
Control 2:	
Control 3:	
Complet Factor:	
Complet Offset:	
DUAL MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
SAMPLE	
Normal Volume:	
Reagent Volume:	
R2 Bottle:	
Normal Volume:	
Reagent Volume:	
Predilution:	
Slope Blank:	
Delay Min Time:	
Linearity Limit:	
Point One Two:	
Incubation Time:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
Substr Depletion:	
Reagent Blank:	
R ABS Deviation:	
Cal Low Limit:	
Cal High Limit:	
Factor:	
MONO MODE	
Name:	
Sample Blank:	
R1 Bottle:	
Normal Volume:	
Reagent Volume:	
SAMPLE	
Normal Volume:	
Reagent Volume:	
Delay Min Time:	
Linearity Limit:	
Predilution:	
Incubation Time:	
Point One Two:	
Low Absorbance:	
High Absorbance:	
R ABS L Limit:	
R ABS H Limit:	
Substr Depletion:	
Reagent Blank:	
R ABS Deviation:	
Cal Low Limit:	
Cal High Limit:	
Factor:	
(*)	
(**)	
(***)	
(#)	

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Name:	
Immunosay:	
Chemistry Type:	
BioType:	
Immuno Chemistry:	
% Sample Volume:	
Wavelength:	
Bichromatic Chemistry:	
Bichromatic Factor:	
K1:	
K2:	
Bio Limit 1:	
Bio Limit 2:	
Deviation Limit:	
Delay Time:	
Blank Type:	
Incubation:	
% Reagent Volume:	
2nd Reagent:	
2nd Reagent Volume:	
A2 Delay:	
Unit:	
Unit Factor:	
Decimal Point:	
RBL Low:	
RBL High:	
Range Low:	
Range High:	
Calibration Factor:	
Reagent Rate:	
Standard Value:	
Normal Low:	
Normal High:	
Slope:	
Intercept:	
C1*10E-6:	
C2*10E-6:	
D1*10E-6:	
Delta R:	
Linear Factor:	
Final Limit:	
Endpoint Limit:	
(*)	
(*)	

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Measure Unit:	Decimals:	Reaction Type:
Volume (µ):	Serum:	Plasma:
	Reagent 1:	Reagent 2:
Abs Range (n Abs):	Min:	Max:
Linearity Limit:		Reagent Blanking:
Contaminating:		Differential:
Filter 1:		Filter 2:
Time (sec):	Mix 1:	Incubation 1:
	Mix 2:	Incubation 2:
Measurement Type:		Factor:
Normal Range:		
Age (yrs):	Min:	Max:
Below 10:		Min:
From 10 to 60:		Max:
Other:		

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