QUIMIFRU – Frutosamina

Colorimétrico

ANVISA: 10159820277

Finalidade

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3. F. San-Gil, G. M. Schier, R. G. Moses, and I. E. Gan, "Improved estimation of fructosamine, as a measure of glycated serum protein, with the Technicon RA-1000 analyzer,," Clin Chem, vol. 31, no. 12, pp. 2005–6, Dec. 1985, Accessed: Feb. 28, 2020. [Online]. Available: <http://www.ncbi.nlm.nih.gov/pubmed/4064290>

4. J. D. Kruse-Jarres, J. Jarausch, P. Lehmann, B. W. Vogt, and P. Rietz, "A new colorimetric method for the determination of fructosamine,," 1989.

5. W. H. O. (WHO), "Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1 Rev.2," 2002. Accessed: Mar. 28, 2019. [Online]. Available: https://apps.who.int/iris/bitstream/handle/10665/65957/WHO_DIL_LAB_99.1_REV.2.pdf?sequence=1&isAllowed=y

6. CLSI, EP 05–A3 Evaluation of Precision of Quantitative Measurement Procedures, vol. 05, no. January. 2014.

7. CLSI, EP17–A2 Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition, Second edi., no. June. 2012.

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8. CLSI, EP07 Interference Testing in Clinical Chemistry, 3rd Edition. 2018.

9. D. S. Young, Effects of drugs on clinical laboratory tests, 5th ed. Washington D.C.: AACC Press, 2000.

10. CLSI, EP 6–A Evaluation of the Linearity of Quantitative Measurement procedures a Statistical Approach. 2003.

11. CLSI, EP25–A Evaluation of Stability of In vitro Diagnostic Reagents; Approved Guideline. 2009.

12. NCCLS-EP9–A2, Method Comparison and Bias Estimation Using Patient samples; Approved Guideline–Second edition NCCLS-EP9–A2, no. 19.

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ANVISA: 10159820277

Purpose

Reagent used for the quantitative determination of fructosamine concentration in human serum or plasma. For in vitro diagnostic use only.

Principle

QUIMIFRU – Fructosamine is based on the reducing capacity of ketoamines formed by the non-enzymatic condensation of glucose and proteins (fructosamines), which reduce tetrazolium salts in an alkaline medium, forming a colored formazan complex.

The rate of color formation is directly proportional to the fructosamine concentration present in the sample.

Methodology

Colorimetric – NBT

Clinical Significance

Fructosamine is the generic term used to describe the ketoamines formed by the non-enzymatic binding of glucose to plasma proteins. The concentration of fructosamine in the blood reflects the average blood glucose level over the previous 2 to 3 weeks.

The determination of fructosamine is important for the short-term monitoring of therapeutic effectiveness in diabetic patients, complementing glycated hemoglobin (HbA1c), which reflects long-term glycemic control.

Fructosamine measurement may also be useful in situations where HbA1c has clinical limitations, such as in patients with altered erythrocyte lifespan. However, significant changes in the concentration or half-life of plasma proteins may affect the results.

Reagents

Reagent: Ready for use. Store at 2–8°C. Carbonate buffer 200 mmol/L, pH 10.35, NBT 0.25 mmol/L, surfactants 2 g/L, and stabilizers.

Calibrator: Ready for use. Store at 2–8°C. Glycated albumin in 0.9 g/L NaCl saline solution.

The reagents are stable until the expiration date printed on the labels when stored refrigerated at 2–8°C, protected from light, and handled in a manner that prevents contamination during use. The on-board stability of the reagent (when stored in the analyzer's refrigerated compartment) is 60 days.

Precautions and Required Care

1. This reagent is for in vitro diagnostic use only.
2. Do not use the reagent and/or calibrator after the expiration date indicated on the packaging.
3. The frequency of assay calibration depends on the instrument used. In addition, we recommend recalibrating the assay and running controls with each new reagent lot.
4. Do not pipette by mouth. Handle with the standard precautions required for laboratory reagents.
5. Samples containing human-derived materials should be handled as potentially infectious using safe laboratory procedures.
6. Adopt standard precautions for handling all laboratory reagents. Infectious or microbiological waste should be treated as potentially hazardous material.
7. Do not ingest or allow contact with skin or eyes. Additional safety information is available in the product Safety Data Sheet (SDS). To obtain the SDS, visit ebram.com.br or contact Ebram customer service.
8. Dispose of in accordance with established laboratory procedures. Environmental hazards: follow all applicable local regulations for safe disposal.

Materials Required but Not Provided

- Water bath or analyzer capable of maintaining a temperature of 37°C and measuring absorbance at 540 nm.
- Pipettes for measuring samples and reagents.
- Distilled/deionized water.
- Analyzer consumables, when required.
- Controls.
- Stopwatch.

Sample

Fresh serum or plasma collected with EDTA or heparin.

Hemolyzed or contaminated samples are not suitable for testing.

Fructosamine is stable for up to 2 weeks at 4–8°C or up to 2 months at –20°C.

Allow samples to reach room temperature before testing. Frozen samples must be completely thawed and thoroughly homogenized before testing. Samples should not be subjected to repeated freeze–thaw cycles.

Patient Preparation

It is recommended that samples be collected in accordance with medical recommendations and/or the reference standard used by the laboratory.

Interferences

Bilirubin above 20 mg/dL and lipemia above 10 g/L may interfere with assay results. Hemoglobin up to 10 g/L does not present significant interference. Other drugs and substances may also interfere with results.

System Parameters

Temperature	37°C
Wavelength	540 nm
Reaction type	Fixed Time
Calibration	Two-point (water and calibrator)
Direction	Increasing
Sample volume	10 µL
Reagent volume	200 µL
Reading 1 (A1)	After 5 minutes from sample addition
Reading 2 (A2)	After 5 minutes from the first reading

Calibration

The method must be calibrated using the fructosamine calibrator provided in the kit, the value of which is traceable to the NCCLS EP9-A2 reference material.

Automated Procedure

Application on automated systems: refer to the equipment operating manual or contact customer support.

Semi-automatic Procedure

Allow reagent, calibrator, and sample to reach room temperature (15–30°C) before testing.

1. Separate 3 test tubes and perform the procedure as shown in the table below:

	Blank	Calibrator	Sample
Reagent (µL)	1000	1000	1000
Deionized water (µL)	50	-	-
Calibrator (µL)	-	50	-
Sample (µL)	-	-	50

2. Mix the tubes and incubate at 37°C.
3. Perform the 1st reading (A1) 5 minutes into the test, zeroing against the reagent blank.
4. Perform the 2nd reading (A2) after an additional 5 minutes of incubation.

Calculations

The fructosamine concentration in the sample is calculated using the following formula:

Sample Rate:

$$\text{Sample Rate} = \frac{[A1] - [A2]}{5}$$

Calibration Rate:

$$\text{Calibration Rate} = \frac{[A1] - [A2]}{5}$$

Sample Calculation:

$$\text{Fructosamine concentration} = \frac{[\text{Sample rate}] - [\text{Calibration rate}]}{[\text{Cal.}]}$$

Linearity

The assay has a measuring range of 50.66 to 800 µmol/L.

Samples with concentrations exceeding the linear range must be diluted with saline solution, and the final result multiplied by the dilution factor.

Quality Control

Each laboratory must maintain an internal quality program that defines objectives, procedures, standards, tolerance limits, and corrective actions. A defined system should also be maintained to monitor analytical variation of the measurement system. The use of commercial controls with manufacturer-assigned target values is recommended.

Reference Values

Serum or plasma: 205–285 µmol/L for the non-diabetic population. In children, values may be approximately 5% lower. These values are provided for guidance only; each laboratory should establish its own reference intervals.

Precision

Repeatability

Repeatability was evaluated using three serum samples with different fructosamine concentrations (Low – 260 µmol/L, Medium – 402.5 µmol/L, and High – 545 µmol/L). Each sample was analyzed 40 times in duplicate on the same day under the same conditions, totaling 80 measurements per level.

n=80	Low	Medium	High
Mean (µmol/L)	260,35	402,23	546,45
C.V. (%)	1,14%	1,09%	1,17%

Reproducibility

Reproducibility was evaluated using three serum samples with different fructosamine concentrations (Low – 260 µmol/L, Medium – 402.5 µmol/L, and High – 545 µmol/L). The samples were analyzed in duplicate over 20 consecutive days, totaling 80 measurements per level

n=80	Low	Medium	High
Mean (µmol/L)	260,34	401,94	545,53
C.V. (%)	1,58%	1,55%	1,67%

Accuracy

The fructosamine assay was tested on the EBIO Plus 200 biochemical analyzer, and the results obtained were compared with a commercial reference method. A total of 47 serum samples, ranging from 110 to 720 µmol/L of fructosamine, were tested in both assays and showed a correlation coefficient of 0.9781. Linear regression analysis resulted in the following equation: y = 1.0085x – 6.8133.

Limit of Detection

The limit of detection (LoD) of QUIMIFRU – Fructosamine was determined using two reagent lots, totaling 120 measurements performed in different matrices containing low fructosamine concentrations. The assay showed a limit of detection of 19.86 µmol/L, with coefficients of variation below 5% across all evaluated samples

Specificity

The analytical specificity of QUIMIFRU – Fructosamine was evaluated using serum samples containing approximately 200 µmol/L and 400 µmol/L of fructosamine, analyzed in the absence and presence of potential interfering substances. The results demonstrated interference from bilirubin at concentrations above 20 mg/dL and lipids above 10 g/L. Hemoglobin up to 10 g/L did not show significant interference in the assay.

Notes

Clinical diagnosis should not be based solely on the results of a single test; therefore, the patient's clinical data as well as the results of other tests must be considered for diagnostic conclusion

Presentation

Reagent: 1 x 50 mL + Cal.: 1 x 1 mL

Bibliographic Reference

1. N. Rifai, A. R. Horvath, and C. T. Wittwer, Tietz textbook of clinical chemistry and molecular diagnostics, 6xt edition. 2018.
2. R. N. Johnson, P. A. Metcalf, and J. R. Baker, "Fructosamine: a new approach to the estimation of serum glycosylprotein. An index of diabetic control," 1982.
3. F. San-Gil, G. M. Schier, R. G. Moses, and I. E. Gan, "Improved estimation of fructosamine, as a measure of glycated serum protein, with the Technicon RA-1000 analyzer," Clin Chem, vol. 31, no. 12, pp. 2005–6, Dec. 1985, Accessed: Feb. 28, 2020. [Online]. Available: <http://www.ncbi.nlm.nih.gov/pubmed/4064290>
4. J. D. Kruse-Jarres, J. Jarausch, P. Lehmann, B. W. Vogt, and P. Rietz, "A new colorimetric method for the determination of fructosamine," 1989.
5. W. H. O. (WHO), "Use of anticoagulants in diagnostic laboratory investigations. Document WHO/DIL/LAB/99.1 Rev.2," 2002. Accessed: Mar. 28, 2019. [Online]. Available: https://apps.who.int/iris/bitstream/handle/10665/65957/WHO_DIL_LAB_99.1_REV.2.pdf?sequence=1&isAllowed=y
6. CLSI, EP 05–A3 Evaluation of Precision of Quantitative Measurement Procedures, vol. 05, no. January, 2014.
7. CLSI, EP17–A2 Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition, Second edi., no. June. 2012.
8. CLSI, EP07 Interference Testing in Clinical Chemistry, 3rd Edition. 2018.
9. D. S. Young, Effects of drugs on clinical laboratory tests, 5th ed. Washington D.C.: AACC Press,

2000.

10. CLSI, EP 6–A Evaluation of the Linearity of Quantitative Measurement procedures a Statistical Approach. 2003.
11. CLSI, EP25–A Evaluation of Stability of In vitro Diagnostic Reagents; Approved Guideline. 2009.
12. NCCLS–EP9–A2, Method Comparison and Bias Estimation Using Patient samples; Approved Guideline–Second edition NCCLS–EP9–A2, no. 19.
13. Arquivos Ebram.

Revision: July 2026

Universal Symbols Used on In Vitro Diagnostic Packaging

	Consult instructions for use		Reagent		Manufactured by
	The content is sufficient for <n> tests		Expiration date (last day of the month)		Lot number
	Temperature limit (store at)		In vitro diagnostic medical device		Catalogue number

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