Quimicalib Calibrador

Calibrador Multiparâmetro

ANVISA: 10159820122

Finalidade

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Referências

1. Tietz Textbook of Clinical Chemistry, edited by Carl A. Burtis and Edward R. Ashwood, W.B Saunders Company, Philadelphia, Second Edition, (1986).

2. Tietz, N. W., Fundamentals of Clinical Chemistry, W. B. Saunders, Philadelphia (1989).

3. Wick, Pinggera W., Lehmann P., Ferritin in Iron Metabolism – Diagnosis of Anemias Springer-Verlag, Second Edition (1995).

4. Morikawa K., Oseko F., Morikawa S., A role for ferritin in hematopoiesis and the immune system. Leuk-Lymphoma 18 (5-6); 429-433, (1995).

5. Domingez T.P., Martinez J.M., Fuentes-Arderiu X., Evaluation of a new measurement procedure for the concentration of ferritin in serum. Eur.J.Clin.Chem.Clin.Biochem., 35 (2); 117-120, (1997).

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ANVISA: 10159820122

Purpose

Multi-parameter biochemistry calibrator, for calibrating clinical chemistry assays in automatic analyzers and manual procedures. "For in vitro diagnostic use only".

Content

The Calibrator is a lyophilized human serum base intended for use as a calibrator in clinical chemistry assays. Check the concentrations of the analytes on the back of this instruction for use. ATTENTION! Concentrations vary per batch.

Presentation

Code 7023 : 5 x 3.0 mL

Code 12023 : 1 x 3.0mL

Traceability

The concentrations of the analytes present in Quimicalib - Calibrator are traceable to standard reference materials which vary according to the analyte. Check the reference material to which each analyte is traceable on the back of this instruction for use.

Preparation instructions

1. Allow the calibrator to reach a temperature of 22 to 28°C.
2. Tap the vial lightly with your fingers to loosen the lyophilized material.
3. Remove the plastic cap and rubber cap from the vial of lyophilized serum.
4. Using a calibrated volumetric pipette, add exactly 3.0 mL of distilled or deionized water to the lyophilized serum.
5. Replace the rubber cap on the vial and leave to stand for 30 minutes.
6. Gently homogenize the contents until the lyophilisate has completely dissolved.
7. Always homogenize gently by swirling before use and then remove an aliquot for testing.

Note: Keep the calibrator open for as short a time as possible, outside the storage temperature or exposed to light.

Storage and Stability

Keep the lyophilized calibrator refrigerated between 2 and 8°C, the unopened product is stable until the expiration date printed on the label.

During handling, the calibrator is subject to contamination of a chemical and microbial nature that can cause a reduction in stability, so to avoid this, it should be handled in accordance with Good Laboratory Practice.

After reconstitution, the analytes in Quimicalib - Calibrator are stable for 7 days if stored between 2 - 8°C in a well-sealed vial and protected from light, except for the Direct Bilirubin, Total Bilirubin which have a stability of 1 day and Alkaline Phosphatase which has a stability of 3 days if stored between 2 - 8°C.

Precautions and Cautions Required

- This calibrator should only be used for "in vitro" diagnostics.
- Do not pipette by mouth. Avoid contact with skin and clothing. In case of contact with eyes, rinse with plenty of water and seek medical attention.

- Chemicalib - Calibrator contains sodium azide as a preservative (<1g/L). This component can react with copper and lead and can become an explosive metal. When disposing of it, add a large quantity of water.
- The temperature of the working environment and the reaction time should be monitored in order to obtain correct results.
- The reagent was tested and the results were negative for the presence of HBsAg and Anti-HCV and anti-HIV antibodies. However, the product should be handled with caution as it is potentially infectious due to the presence of human blood derivatives.
- To dispose of the product, we suggest applying local, state or federal environmental protection regulations.

Limitations

- This product should not be used after the expiration date
- Use a calibrated pipette for reconstitution, as a volume other than 3.0 mL may make the calibrator too dilute or concentrated and may alter the pre-established concentrations of the analytes present.
- Some variables can interfere with the results obtained with the calibrator, such as: water quality, contamination of the pipette or glassware, incorrect homogenization, errors in reconstitution, technical faults with the analyzer or reagents used. We suggest following the requirements for good laboratory practice.
- The concentration of the analytes in Chemicalib - Calibrator were established using the Ebram procedure and reagents as described in the instructions for use. The accuracy of the calibration may not be adequate when reagents from other manufacturers are used.

Note

- Proper cleaning and drying of the material used are fundamental factors for the stability of the calibrator and obtaining correct results.
- The water used in the laboratory must be of the appropriate quality for each application. Thus, for preparing reagents and for final rinsing of glassware, the water must have a resistivity greater than or equal to 1 megohm.cm or a conductivity less than or equal to 1 microsiemens/cm and a silicate concentration < 0.1 mg/L.

Quality Assurance

This product is guaranteed by Ebram Produtos Laboratoriais Ltda if stored at the recommended temperature, used within the expiration date and following the recommendations on the label and in this instruction for use.

References

1. Tietz Textbook of Clinical Chemistry, edited by Carl A. Burtis and Edward R. Ashwood, W.B Saunders Company, Philadelphia, Second Edition, (1986).
2. Tietz, N. W., Fundamentals of Clinical Chemistry, W. B. Saunders, Philadelphia (1989).

3. Wick, Pinggera W., Lehmann P., Ferritin in Iron Metabolism – Diagnosis of Anemias Springer-Verlag, Second Edition (1995).
4. Morikawa K., Oseko F., Morikawa S., A role for ferritin in hematopoiesis and the immune system. Leuk-Lymphoma 18 (5-6); 429-433, (1995).
5. Domingez T.P., Martinez J.M., Fuentes-Arderiu X., Evaluation of a new measurement procedure for the concentration of ferritin in serum. Eur.J.Clin.Chem.Clin.Biochem., 35 (2); 117-120, (1997).

Revision: March 2026.

Universal symbols used on packaging in vitro diagnostics			
	See instructions for use	 Reagent	 Manufactured by
	The contents are sufficient for <n> tests	 Expiry date (last day of the month)	 Batch number
	Temperature limit (store at)	 Product for in vitro diagnostics	 Catalogue number
	Use veterinarian		

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Responsible Technician
 Dr. Nadjara Novaes Longen – CRF-SP – 37.451

SAC Ebram
 For further information, please contact
 SAC 0800 500 2424 or ☎ 11 2574 7110
 sac@ebram.com | ebram.com.br