# Quimicontrol

## Soro Controle Anormal

ANVISA: 10159820155

### Finalidade

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ANVISA: 10159820155

### Purpose

The Abnormal Control Serum is intended for use in internal quality control in clinical laboratories and is suitable for evaluating the accuracy, precision and calibration deviations of equipment used for quantitative dosages in clinical chemistry.

For "in vitro" diagnostic use only.

### Presentation

Code: 7031 - Abnormal Control Serum: 10 X 5.0ml.

Code: 12031 - Abnormal Control Serum: 3 X 5.0ml.

Code: 12051 - Abnormal Control Serum: 1 X 5,0 mL

### Composition

Lyophilized serum with various analytes that are adjusted by the addition of animal tissue extracts and other non-protein constituents such as purified chemical materials.

### Characteristics

Controls are used to evaluate the reliability of an analytical process, particularly with regard to precision, accuracy and also as a means of monitoring analytical performance day by day, month by month, etc. GLP (Good Clinical Laboratory Practice) recommends that laboratories use at least 2 controls with different concentration levels, so that the information is valid for verifying that the desired control levels are maintained, thus minimizing unnecessary repetitions caused by false rejection, which is commonly associated with the use of a single level of control serum.

Through appropriate management decisions, the information generated through the use of controls can be used to verify adherence to the necessary performance requirements.

### Cautions

The serum used in this product has been tested and the results obtained were negative for HBsAg and anti-HIV, but it is still recommended that this product be handled with the same precautions used with patient samples.

### Preparation

\*Allow the control to reach room temperature before reconstituting.

1. Remove the rubber cap from the vial to be used.
2. Using a volumetric pipette, add exactly 5mL of distilled or deionized water to the vial.
3. Replace the cap on the vial and leave to stand for 30 minutes.
4. Shake the contents gently until dissolved. Avoid foaming.
5. Take out the amount required for use.
6. Cap immediately and store away from light at 2-8°C.

The Control Serum is to be used in accordance with the guidelines accompanying the test procedure and in the same way as the patient samples.

### Stability and Conservation

**Unreconstituted material:** The Control Serum should be stored at 2 - 8°C and is stable until the expiration date indicated on the vial labels and packaging.

**Reconstituted material:** constituents are stable for up to 7 days at 2 - 8°C or for up to 28 days at -20°C, when stored in a well-sealed vial and protected from light, except for the analytes Alkaline Phosphatase and Bilirubin, whose stability is 2 days, when stored at 2 - 8°C with the vial well sealed.

The control should be protected from light and kept under refrigeration when not in use.

To avoid freezing and thawing, we suggest separating the control serum into aliquots of 0.5mL or 1mL in hermetically sealed vials.

Before use, the product kept refrigerated or frozen should be properly homogenized.

### Recommendation

- We recommend using EBRAM Normal Control Serum and EBRAM Abnormal Control Serum together in order to obtain greater precision in the Quality Control carried out.
- The concentration ranges described in this instruction for use were established in an analytical system where the reagents and equipment used are manufactured by Ebram. We suggest that the laboratory establishes its own average value for the analytes, taking into account the laboratory's analytical variability (equipment, water, inaccuracies in measurement procedures, etc.).
- The laboratory's Technical Manager or Quality Control Manager must evaluate the daily dosages of this control serum and make the appropriate decisions for the release of patient results.
- The water used in the laboratory must be of the appropriate quality for each application. Thus, to prepare reagents and use in measurements, it must have resistivity  $\geq 1$  mega ohm or conductivity  $\leq 1$  microsiemens and silicate concentration  $< 0.1$  mg/L (purified water). For rinsing the glassware, the water can be purified, with resistivity  $\geq 0.1$  megaohms or conductivity  $\leq 10$  microsiemens. Use purified water for the final rinse.

### Limitations

- Failure to obtain expected values may indicate inadequate reconstitutions, the need for analyzer maintenance or errors in technique associated with the test procedure.
- Improper storage or handling of control serum can also affect the results, do not use vials with visible evidence of microbial growth.
- The ranges were determined after statistical treatment of the results obtained with the EBRAM reagents on specific equipment, so differences in the analytical system may result in different results.
- Thorough washing and drying of the glassware, checking the stability of the reagents, correct pipetting, observing the reaction time and temperature are extremely important for obtaining precise and accurate results.

Revision: March 2026.

#### Universal symbols used on packaging in vitro diagnostics

|                                                                                                                                 |                                                                                                                           |                                                                                                        |
|---------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
|  See instructions for use                    |  Reagent                             |  Manufactured by  |
|  The contents are sufficient for $<n>$ tests |  Expiry date (last day of the month) |  Batch number     |
|  Temperature limit (store at)                |  Product for in vitro diagnostics    |  Catalogue number |
|  Use veterinarian                            |                                                                                                                           |                                                                                                        |

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