Quimicontrol

Soro Controle Normal

ANVISA: 10159820155

Finalidade

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ANVISA: 10159820155

Purpose

Normal Control Serum is intended for use in internal quality control in clinical laboratories, and is suitable for evaluating the accuracy, precision and calibration deviations of the equipment used for quantitative dosages in clinical chemistry. For "in vitro" diagnostic use only

Presentation

Code: 7024 - Normal Control Serum: 10 x 5.0mL

Code: 12024 - Normal Control Serum: 3 x 5.0mL

Code: 12044 - Normal Control Serum: 1 x 5.0mL

Composition

Lyophilized serum with various analytes that are adjusted by the addition of animal tissue extracts and other non-protein constituents such as purified chemical materials.

Characteristics

Controls are used to evaluate the reliability of an analytical process, particularly with regard to precision, accuracy and also as a means of monitoring analytical performance day by day, month by month, etc. GLP (Good Clinical Laboratory Practice) recommends that laboratories use at least 2 controls with different concentration levels, so that the information is valid for verifying that the desired control levels are being maintained, thus minimizing unnecessary repetitions caused by false rejection, which is commonly associated with the use of a single level of control serum.

Through appropriate management decisions, the information generated through the use of controls can be used to verify adherence to the necessary performance requirements.

Cautions

The serum used in this product has been tested and the results obtained were negative for HBsAg and anti-HIV. As no method can provide complete assurance that a product derived from human blood will not transmit an infectious agent, it is recommended that this product be handled with the same precautions used with patient samples.

Preparation

*Allow the control to reach room temperature before reconstituting.

1. Remove the rubber cap from the vial to be used.
2. Using a volumetric pipette, add exactly 5mL of distilled or deionized water to the vial.
3. Replace the cap on the vial and leave to stand for 30 minutes.
4. Shake the contents gently until dissolved. Avoid the formation of foam.
5. Take out the quantity required for use.
6. Cap immediately and store away from light at 2-8°C.

The Control Serum is to be used in accordance with the guidelines accompanying the test procedure and in the same way as the patient samples.

Stability and Conservation

Unreconstituted material: Control Serum should be stored at 2 - 8°C and is stable until the expiry date indicated on the vial labels and packaging.

Reconstituted material: constituents are stable for up to 7 days when stored in a well-sealed vial at 2 - 8°C and protected from light, except for Alkaline Phosphatase and bilirubin, which are stable for 2 days. The control should be protected from light and kept refrigerated when not in use.

The analytes are stable for 28 days when stored at - 20°C in tightly sealed vials, with the exception of Bilirubin and Alkaline Phosphatase. To avoid freezing and thawing, it is suggested to separate the control serum into 0.5mL or 1mL aliquots in hermetically sealed vials. Before use, the refrigerated or frozen product should be properly homogenized.

Recommendation

- We recommend using EBRAM Normal Control Serum and EBRAM Abnormal Control Serum together to obtain greater precision in the Quality Control performed.
- The concentration ranges described in this instruction for use were established in an analytical system where the reagents and equipment used are manufactured by Ebram. We suggest that the laboratory establishes its own average value for the analytes, taking into account the laboratory's analytical variability (equipment, water, inaccuracies in measurement procedures, etc.).
- The laboratory's Technical Manager or Quality Control Manager must evaluate the daily dosages of this control serum and make the appropriate decisions for the release of patient results.
- The water used in the laboratory must be of the appropriate quality for each application. Thus, to prepare reagents and use in measurements, it must have resistivity ≥ 1 mega ohm or conductivity ≤ 1 microsiemens and silicate concentration < 0.1 mg/L (purified water). For rinsing the glassware, the water can be purified, with resistivity ≥ 0.1 megaohms or conductivity ≤ 10 microsiemens. Use purified water for the final rinse.

Limitations

- Failure to obtain expected values may indicate inadequate reconstitutions, the need for analyzer maintenance or errors in technique associated with the test procedure.
- Improper storage or handling of control serum can also affect the results, do not use vials with visible evidence of microbial growth.
- The ranges were determined after statistical treatment of the results obtained with the EBRAM reagents in specific equipment, so differences in the analytical system may result in different results.
- Thorough washing and drying of the glassware, checking the stability of the reagents, correct pipetting, observing the reaction time and temperature are extremely important for obtaining precise and accurate results.

Revision: October 2025.

Universal symbols used on packaging in vitro diagnostics

 See instructions for use	 Reagent	 Manufactured by	 Contents are sufficient for $<n>$ tests
 Temperature limit (store at)	 Expiry date (last day of the month)	 Batch number	
 Use veterinary	 Product for in vitro diagnostics	 Catalog number	

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