QUIMIURIC – ÁCIDO ÚRICO

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	AURE	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	ÁCIDO ÚRICO
UNIT	mg/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	5	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 25	UPPER	30
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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CALIBRATION:

REQUEST RULE	LINEAR	DIFFERENCE	0
REPLICATES	2	BLANK RESPONSE	0 – 0
INTERVAL	0	COEFF DIFFERENCE	0
SENSITIVITY	0	NON LINEAR SD	0
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	ALB	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	ALBUMINA
UNIT	g/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	630	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 10	UPPER	6
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	AMI	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	AMILASE
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	405	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	5	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	3 / 15	UPPER	1317
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	BT	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	BILI TOTAL
UNIT	mg/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	12	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	120	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	10	LOWER	0
REACTION	-2 / 25	UPPER	20
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	CAL	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	CÁLCIO
UNIT	mg/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	630	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	4	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 10	UPPER	15
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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CALIBRATION:

REQUEST RULE	ONE POINT LINEAR	DIFFERENCE	0
REPLICATES	2	BLANK RESPONSE	0-0
INTERVAL	0	COEFF DIFFERENCE	0
SENSITIVITY	0	NON LINEAR SD	0
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	CKMB	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	CKMB
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	405	UPPER	-
SAMPLE VOL	10	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	25 / 50	UPPER	1000
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	CLO	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	CLORO
UNIT	mmol/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	70
REACTION	0 / 8	UPPER	150
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	COL	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	COLESTEROL
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 25	UPPER	750
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	CRE	LOWER	0
METHOD	FIXED TIME	UPPER	0
DIRECTION	ASCENDING	FULL NAME	CREATININA
UNIT	mg/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	20	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	3 / 13	UPPER	20
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	FERRO	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	FERRO
UNIT	µg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	45	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	50	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	2	LOWER	0
REACTION	-2 / 25	UPPER	1000
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	ONE POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	FAL	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	FAL
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	405	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	4	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	5 / 11	UPPER	1500
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	ONE POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	FÓSFORO	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	FÓSFORO
UNIT	mg/dL	PRINT N°	
DECIMALS	1	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	4	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 25	UPPER	15
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	ONE POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	GLIOX	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	GLICOSE
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 20	UPPER	400
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	HDL	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	HDL COLESTEROL
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	578	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	225	LOWER	-
R2 VOL	75	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	25	LOWER	0
REACTION	-2 / 25	UPPER	180
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	DHL	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	DESCENDING	FULL NAME	DHL
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	5	MIX. R BLANK	
R1 VOL	160	LOWER	-
R2 VOL	40	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	2	LOWER	0
REACTION	5 / 15	UPPER	1250
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	DHL	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	DESCENDING	FULL NAME	DHL
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	700	UPPER	-
SAMPLE VOL	5	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	5 / 15	UPPER	1250
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	LACTE	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	LACTATO
UNIT	mg/dL	PRINT N°	
DECIMALS	1	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0.3
REACTION	0 / 25	UPPER	90
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	MAGNÉSIO	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	MAGNÉSIO
UNIT	mg/dL	PRINT N°	
DECIMALS	2	R1 BLANK	
PRIM WAVE	578	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	4	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	-1 / 8	UPPER	6.1
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	ONE POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	PT	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	PROTEÍNA
UNIT	g/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 23	UPPER	12
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	PTUE	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	PROTEINURIA
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	600	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	6	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	8
REACTION	1 / 19	UPPER	400
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	TGO	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	DESCENDING	FULL NAME	TGO
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	20	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	3 / 13	UPPER	400
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	TGP	LOWER	0
METHOD	KINECT	UPPER	0
DIRECTION	DESCENDING	FULL NAME	TGP
UNIT	U/L	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	20	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	3 / 13	UPPER	400
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	TRI	LOWER	0
METHOD	ENDPOINT	UPPER	0
DIRECTION	ASCENDING	FULL NAME	TRIGLICÉRIDES
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	510	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	0 / 20	UPPER	800
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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PARAMETERS:

N°	(*)	RESPONSE	
TEST	URE	LOWER	0
METHOD	FIXED TIME	UPPER	0
DIRECTION	DESCENDING	FULL NAME	UREIA
UNIT	mg/dL	PRINT N°	
DECIMALS	0	R1 BLANK	
PRIM WAVE	340	LOWER	-
SEC WAVE	-	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	300	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	0
REACTION	3 / 13	UPPER	200
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0-0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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