## QUIMIURIC – ÁCIDO ÚRICO

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**PARAMETERS**

TEST NAME:	AUR	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Ácido Úrico	SAMPLE VOLUME:	5 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 18	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	ADA	R1:	180 µL
TEST N°:		R2:	90 µL
FULL NAME:	Adenosina Deaminase	SAMPLE VOLUME:	6 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	FIXED TIME	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	700	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	18 30	FACTOR:	---
INCUBA TIME:	9	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTERGER		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
1		SERUM	1 - 100 YEAR	0.00	15.0
2		LIQ PLEURAL	1 - 100 YEAR	0.00	30.0
3		LCR	1 - 100 YEAR	0.00	9.0

### CALIBRATION RULE

RULE:	ONE-POINT LINEAR	NAME	CONCENTRATION
SENSITIVITY:		CAL ADA	(#)
REPLICATES:	2		
INTERVAL:			
DIFFERENCE LIMIT:			
SD:			
BLANK RESPONSE:			
ERROR LIMIT:			
COEFFICIENT:			

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### PARAMETERS

TEST NAME:	ALB	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Albumina	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	630	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 4	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	g/dL		
PRECISION:	0.01		

### REFERENCE

N° (*)	GENDER (*)	SAMPLE TYPE (*)	AGE (*)	LOW LIM (*)	UP LIM (*)
-----------	---------------	--------------------	------------	----------------	---------------

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	AMI	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Amilase	SAMPLE VOLUME:	5 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINECT	MIXED REAG BLK:	---
PRI. WAVE:	405	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	BD	R1:	200 µL
TEST N°:		R2:	50 µL
FULL NAME:	Bili Direta	SAMPLE VOLUME:	25 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 13	FACTOR:	---
INCUBA TIME:	4	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	BT	R1:	200 µL
TEST N°:		R2:	120 µL
FULL NAME:	Bili Total	SAMPLE VOLUME:	12 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 19	FACTOR:	---
INCUBA TIME:	4	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	CAL	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Cálcio	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	630	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 8	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	CKMB	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	CKMB	SAMPLE VOLUME:	8 µL
STANDARD N°:		R1 BLANK:	---
ANALY. TYPE:	KINECT	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	22 34	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	

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**PARAMETERS**

TEST NAME:	CKNAC	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	CKNAC	SAMPLE VOLUME:	10 µL
STANDARD N°:		R1 BLANK:	---
ANALY. TYPE:	KINECT	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 12	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	CLO	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	COLORO	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 8	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mmol/L		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	COL	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Colesterol	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 25	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	CHE	R1:	200 µL
TEST N°:		R2:	40 µL
FULL NAME:	CHE	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINETIC	MIXED REAG BLK:	---
PRI. WAVE:	405	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	DESCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	24 36	FACTOR:	102000
INCUBA TIME:	19	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	CRE	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Creatinina	SAMPLE VOLUME:	20 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	FIXED TIME	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	3 9	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	FER	R1:	160 µL
TEST N°:		R2:	40 µL
FULL NAME:	Ferro	SAMPLE VOLUME:	40 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 32	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	ug/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	FAL	R1:	250 µL
TEST N°:		R2:	0 µL
FULL NAME:	Fosfatase Alcalina	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINECT	MIXED REAG BLK:	---
PRI. WAVE:	405	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	FOS	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Fósforo	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 18	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	GGT	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Gama GT	SAMPLE VOLUME:	8 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINECT	MIXED REAG BLK:	---
PRI. WAVE:	405	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	GLI	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Glicose	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 20	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	HDL	R1:	300 µL
TEST N°:		R2:	100 µL
FULL NAME:	HDL Colesterol	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	578	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 19	FACTOR:	---
INCUBA TIME:	18	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	DHL	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	DHL	SAMPLE VOLUME:	5 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINETIC	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	DESCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	LACTE	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Lactato	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	0 19	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N° (*)	GENDER (*)	SAMPLE TYPE (*)	AGE (*)	LOW LIM (*)	UP LIM (*)
-----------	---------------	--------------------	------------	----------------	---------------

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	LIPAE	R1:	200 µL
TEST N°:		R2:	50 µL
FULL NAME:	Lipase	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINETIC	MIXED REAG BLK:	---
PRI. WAVE:	578	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	0 19	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	1		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	MAG	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Magnésio	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	578	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 8	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	PROT	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Proteína	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	546	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 23	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	g/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	PTUE	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	Proteinúria	SAMPLE VOLUME:	4 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	600	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	1 19	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	TGO	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	TGO	SAMPLE VOLUME:	20 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINETIC	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	DESCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	TGP	R1:	200 µL
TEST N°:		R2:	0 µL
FULL NAME:	TGP	SAMPLE VOLUME:	20 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	KINETIC	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	DESCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	5 11	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	U/L		
PRECISION:	INTEGER		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	ONE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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### PARAMETERS

TEST NAME:	TRI	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Triglicérides	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	ENDPOINT	MIXED REAG BLK:	---
PRI. WAVE:	510	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	ASCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	-1 20	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

### REFERENCE

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

### CALIBRATION RULE

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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**PARAMETERS**

TEST NAME:	URE	R1:	300 µL
TEST N°:		R2:	0 µL
FULL NAME:	Ureia	SAMPLE VOLUME:	3 µL
STANDARD N°:	(*)	R1 BLANK:	---
ANALY. TYPE:	FIXED TIME	MIXED REAG BLK:	---
PRI. WAVE:	340	CONCENTRATION:	(@) (@)
SEC. WAVE:	---	LINEARITY LIMIT:	---
TREND:	DESCENDING	SUBSTRATE LIMIT:	---
REAC. TIME:	3 9	FACTOR:	---
INCUBA TIME:	---	PROZONE CHECK:	---
UNIT:	mg/dL		
PRECISION:	0.01		

**REFERENCE**

N°	GENDER	SAMPLE TYPE	AGE	LOW LIM	UP LIM
(*)	(*)	(*)	(*)	(*)	(*)

**CALIBRATION RULE**

RULE:	SINGLE POINT LINEAR
SENSITIVITY:	
REPLICATES:	2
INTERVAL:	
DIFFERENCE LIMIT:	
SD:	0
BLANK RESPONSE:	
ERROR LIMIT:	
COEFFICIENT:	0
CONCENTRATION:	(#)

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