TURB HBA1C – HEMOGLOBINA GLICADA

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PARAMETERS

| | | | |
|--------------|------------|------------------|---------|
| TEST NAME: | HBA1C | R1: | 180 µL |
| TEST N°: | | R2: | 60 µL |
| FULL NAME: | HBA1C | SAMPLE VOLUME: | 12 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | Fixed Time | MIXED REAG BLK: | --- |
| PRI. WAVE: | 630 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | 2 20 | FACTOR: | --- |
| INCUBA TIME: | 18 | PROZONE CHECK: | --- |
| UNIT: | % | | |
| PRECISION: | 0.01 | | |

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PARAMETERS

| | | | |
|--------------|----------------|------------------|---------|
| TEST NAME: | PCTE | R1: | 180 µL |
| TEST N°: | | R2: | 60 µL |
| FULL NAME: | Procalcitonina | SAMPLE VOLUME: | 20 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | Fixed Time | MIXED REAG BLK: | --- |
| PRI. WAVE: | 578 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | 6 19 | FACTOR: | --- |
| INCUBA TIME: | 17 | PROZONE CHECK: | --- |
| UNIT: | ng/mL | | |
| PRECISION: | 0.01 | | |

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PARAMETERS

| | | | |
|--------------|-----------|------------------|---------|
| TEST NAME: | ASO | R1: | 240 µL |
| TEST N°: | | R2: | 0 µL |
| FULL NAME: | ASO | SAMPLE VOLUME: | 3 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | ENDPOINT | MIXED REAG BLK: | --- |
| PRI. WAVE: | 546 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | -1 10 | FACTOR: | --- |
| INCUBA TIME: | --- | PROZONE CHECK: | --- |
| UNIT: | UI/mL | | |
| PRECISION: | 0.01 | | |

REFERENCE

| N° | GENDER | SAMPLE TYPE | AGE | LOW LIM | UP LIM |
|-----|--------|-------------|-----|---------|--------|
| (*) | (*) | (*) | (*) | (*) | (*) |

CALIBRATION RULE

| | |
|-------------------|---------------------|
| RULE: | SINGLE POINT LINEAR |
| SENSITIVITY: | |
| REPLICATES: | 2 |
| INTERVAL: | |
| DIFFERENCE LIMIT: | |
| SD: | 0 |
| BLANK RESPONSE: | |
| ERROR LIMIT: | |
| COEFFICIENT: | |
| CONCENTRATION: | (#) |

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PARAMETERS

| | | | |
|--------------|-----------|------------------|---------|
| TEST NAME: | PCR Ultra | R1: | 225 µL |
| TEST N°: | | R2: | 0 µL |
| FULL NAME: | PCR Ultra | SAMPLE VOLUME: | 3 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | ENDPOINT | MIXED REAG BLK: | --- |
| PRI. WAVE: | 546 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | -1 18 | FACTOR: | --- |
| INCUBA TIME: | --- | PROZONE CHECK: | --- |
| UNIT: | mg/L | | |
| PRECISION: | 0.01 | | |

REFERENCE

| N° | GENDER | SAMPLE TYPE | AGE | LOW LIM | UP LIM |
|-----|--------|-------------|-----|---------|--------|
| (*) | (*) | (*) | (*) | (*) | (*) |

CALIBRATION RULE

| | |
|-------------------|---------------------|
| RULE: | SINGLE POINT LINEAR |
| SENSITIVITY: | |
| REPLICATES: | 2 |
| INTERVAL: | |
| DIFFERENCE LIMIT: | |
| SD: | 0 |
| BLANK RESPONSE: | |
| ERROR LIMIT: | |
| COEFFICIENT: | |
| CONCENTRATION: | (#) |

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PARAMETERS

| | | | |
|--------------|-----------|------------------|---------|
| TEST NAME: | PCR | R1: | 300 µL |
| TEST N°: | | R2: | 0 µL |
| FULL NAME: | PCR | SAMPLE VOLUME: | 3 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | ENDPOINT | MIXED REAG BLK: | --- |
| PRI. WAVE: | 546 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | -1 10 | FACTOR: | --- |
| INCUBA TIME: | --- | PROZONE CHECK: | --- |
| UNIT: | mg/L | | |
| PRECISION: | 0.01 | | |

REFERENCE

| N° | GENDER | SAMPLE TYPE | AGE | LOW LIM | UP LIM |
|-----|--------|-------------|-----|---------|--------|
| (*) | (*) | (*) | (*) | (*) | (*) |

CALIBRATION RULE

| | |
|-------------------|---------------------|
| RULE: | SINGLE POINT LINEAR |
| SENSITIVITY: | |
| REPLICATES: | 2 |
| INTERVAL: | |
| DIFFERENCE LIMIT: | |
| SD: | 0 |
| BLANK RESPONSE: | |
| ERROR LIMIT: | |
| COEFFICIENT: | |
| CONCENTRATION: | (#) |

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PARAMETERS

| | | | |
|--------------|------------|------------------|---------|
| TEST NAME: | VITDE | R1: | 180 µL |
| TEST N°: | | R2: | 45 µL |
| FULL NAME: | Vitamina D | SAMPLE VOLUME: | 3.4 µL |
| STANDARD N°: | (*) | R1 BLANK: | --- |
| ANALY. TYPE: | Fixed Time | MIXED REAG BLK: | --- |
| PRI. WAVE: | 700 | CONCENTRATION: | (@) (@) |
| SEC. WAVE: | --- | LINEARITY LIMIT: | --- |
| TREND: | ASCENDING | SUBSTRATE LIMIT: | --- |
| REAC. TIME: | 2 16 | FACTOR: | --- |
| INCUBA TIME: | 21 | PROZONE CHECK: | --- |
| UNIT: | ng/mL | | |
| PRECISION: | 0.1 | | |

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