## TURB HBA1C – HEMOGLOBINA GLICADA

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### CALIBRATION:

REQUEST RULE	SPLINE	DIFFERENCE	0
REPLICATES		BLANK RESPONSE	0 – 0
INTERVAL	0	COEFF DIFFERENCE	0
SENSITIVITY	0	NON LINEAR SD	0
CORRELATION	0		

CALIBRATION	
Calibrador 1	(#)
Calibrador 2	(#)
Calibrador 3	(#)
Calibrador 4	(#)

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<b>CALIBRATION:</b>		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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**PARAMETERS:**

Nº	(*)	RESPONSE	
TEST	PCRUS	LOWER	0
METHOD	FIXED TIME	UPPER	0
DIRECTION	ASCENDING	FULL NAME	PCRUS
UNIT	mg/L	PRINT Nº	
DECIMALS	0	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	3	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	(@)
REACTION	2 / 34	UPPER	(@)
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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**PARAMETERS:**

Nº	(*)	RESPONSE	
TEST	PCR	LOWER	0
METHOD	FIXED TIME	UPPER	0
DIRECTION	ASCENDING	FULL NAME	PCR
UNIT	mg/L	PRINT Nº	
DECIMALS	2	R1 BLANK	
PRIM WAVE	546	LOWER	-
SEC WAVE	670	UPPER	-
SAMPLE VOL	2	MIX. R BLANK	
R1 VOL	200	LOWER	-
R2 VOL	-	UPPER	-
LINE. LIMIT	-	LINEARITY	
INCUBATION	-	LOWER	(@)
REACTION	2 / 34	UPPER	(@)
ANTIGEN	-	FACTOR	0
SUBSTRATE	-	SPEED	0
CALIBRATION:		DIFFERENCE	0
REQUEST RULE	TWO POINT LINEAR	BLANK RESPONSE	0 – 0
REPLICATES	2	COEFF DIFFERENCE	0
INTERVAL	0	NON LINEAR SD	0
SENSITIVITY	0		
CORRELATION	0		

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**CALIBRATION:**

REQUEST RULE	SPLINE	DIFFERENCE	0
REPLICATES	2	BLANK RESPONSE	0 – 0
INTERVAL	0	COEFF DIFFERENCE	0
SENSITIVITY	0	NON LINEAR SD	0
CORRELATION	0		

CALIBRATION	
VITD CAL 1	(#)
VITD CAL 2	(#)
VITD CAL 3	(#)
VITD CAL 4	(#)
VITD CAL 5	(#)

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