TURB HBA1C – HEMOGLOBINA GLICADA

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Y=ax+b Factor: Direction: Positive

a 1.00
b 0.00

Test of Calibration: 1 STD Number: 4
STD Position and Value: Calibration rules: Spline
#

REAGENT, SAMPLE AND TIME SETUP

R1 Setup:
Volume (uL): 180.0 Position: * Incubation Time (s): 300

R2 Setup:
Volume (uL): 60.0 Position: * Incubation Time (s): 40

Sample Vol. (uL): 12 Test Point (Test Time = Test Point x Test Cycle): 10

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Name: ASOE Byname: ASOE Item Code: *

METHODOLOGY

Assay: End Point Blank Setup: Before adding sample Wave 1: 546 Wave 2: 700
Decimal: 1 Unit: UI/mL Priority: *
☐ Clean Before Test
☐ Clean After Test

Y=ax+b Factor:
a 1.00 BeginPoint (n): 1
b 0.00 EndPoint (m): 2

Test of Calibration: 1 STD Number: 1
STD Position and Value: * Calibration rules: 1-point linear

REAGENT, SAMPLE AND TIME SETUP

R1 Setup:
Volume (uL): 240 Position: * Incubation Time (s): 410
R2 Setup:
Volume (uL): Position: Incubation Time (s):

Sample Vol. (uL): 3.0 Test Point (Test Time = Test Point x Test Cycle): 4

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Name: PCRUSE **Byname:** PCRUSE **Item Code:** *

METHODOLOGY

Assay:	Blank Setup:	Wave 1:	Wave 2:	
End Point	Before adding sample	546	700	
Decimal:	0	Unit:	mg/L	Priority: *
				☐ Clean Before Test
				☐ Clean After Test

Y=ax+b		Factor:		
a	1.00			BeginPoint (n): 1
b	0.00			EndPoint (m): 2

Test of Calibration:	1	STD Number:	1
STD Position and Value:		Calibration rules:	1-point linear
*			

REAGENT, SAMPLE AND TIME SETUP

R1 Setup:		
Volume (uL):	Position:	Incubation Time (s):
240	*	410
R2 Setup:		
Volume (uL):	Position:	Incubation Time (s):

Sample Vol. (uL):	Test Point (Test Time = Test Point x Test Cycle):
3.6	4

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Name: PCRE Byname: PCRE Item Code: *

METHODOLOGY

Assay: Blank Setup: Wave 1: Wave 2:
End Point Before adding sample 546 700
Decimal: 1 Unit: mg/L Priority: *
☐ Clean Before Test
☐ Clean After Test

Y=ax+b Factor:
a 1.00 BeginPoint (n): 1
b 0.00 EndPoint (m): 2

Test of Calibration: 1 STD Number: 1
STD Position and Value: Calibration rules: 1-point linear
*

REAGENT, SAMPLE AND TIME SETUP

R1 Setup:
Volume (uL): Position: Incubation Time (s):
300 * 410

R2 Setup:
Volume (uL): Position: Incubation Time (s):

Sample Vol. (uL): Test Point (Test Time = Test Point x Test Cycle):
3.0 4

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